


# Snoezelen® Balls-in-Bubble Tubes

Enjoy the quiet action of bubbles rising, balls rising and falling, and colored lights changing! This Balls-in-Bubble Tube features a column of rising air bubbles that carry brightly colored balls to the top of the inner tube, gently pushing them over the edge to drift back to the base.



## FEATURES & BENEFITS

- Bubble tubes are used as a focal tool for sensory regulatory activities or used interactively with encouraging visual tracking, color recognition, physical movement (reaching, looking up) and hand-eye coordination and communication
- The balls move slightly slower in this specific bubble tube series, making tracking easier for those with a limited or slower eye tracking
- 8 colors provide a calming, relaxing ambiance
- Bubbles of air propel multi-colored balls to the top of the tube
- Attracts attention, concentration and engagement
- Smooth surface and cooler temperature of the tube provides tactile contrast
- Easy to clean, maintain, operate and assemble
- Made of solid, durable clear acrylic
- Sensory Magic® compatible

## PRODUCT RECOMMENDATIONS

- We recommend using distilled water and adding BCB Water Treatment Fluid (item # 2123508) to the bubble tube monthly or immediately after a water change to help prevent bacteria, smell, and algae
- Top off bubble tube every 3 months, or as needed, with distilled water
- Water level must be kept at least 2" above inner tube to operate properly
- Bubble tube collar is recommended to secure the bubble tube to a wall
- Use in well-lit spaces or in darkened environments where their effects are most striking

## MODELS

MODEL	ITEM #	PRODUCT DIMENSIONS	BASE
70 Inch Ball-in-Bubble Tube	2123704	70"H x 8" Diam	20" Diam
80 Inch Ball-in-Bubble Tube	2123703	80"H x 8" Diam	20" Diam
BUBBLE TUBE BASES			
Square	5011887	40"W X 40"D x 15"H	Additional bubble tube hole variations online
Square	5011888	48"W X 48"D x 15"H	Additional bubble tube hole variations online
Square	5011891	57"W X 57"D x 15"H	Additional bubble tube hole variations online
Triangle	5011895	40"W X 40"D x 15"H	One hole option
Triangle	5011898	48"W X 48"D x 15"H	One hole option
Triangle	5011896	57"W X 57"D x 15"H	Additional bubble tube hole variations online
Quadrant	5011900	40"W X 40"D x 15"H	One hole option
Quadrant	5011899	48"W X 48"D x 15"H	Additional bubble tube hole variations online
Quadrant	5011903	57"W X 57"D x 15"H	Additional bubble tube hole variations online
Half Moon	5011892	48"W X 24"D x 15"H	One hole option
Half Moon	5011894	57"W X 26½"D x 15"H	Additional bubble tube hole variations online

## OPTIONAL ACCESSORIES

MODEL	ITEM #	PRODUCT DIMENSIONS
Sensory Magic®	2123708	
Bubble Tube Collar	2122002	8" Diam
Bubble Tube Filling Pump	2140489	
BCB Water Treatment Fluid	2123508	
Wi-Fi Wireless Controller	2120147	14" Diam x 4¾"H
Interactive Controller	2119933	20" Diam x 4"H
WiFi Interactive Stepping Stones	2120598	57"W x 29"D x 33"H

## CLEANING INSTRUCTIONS

Bubble tubes should ONLY be cleaned with a mild solution of soapy water as described in the instruction sheet. Never use solvent based cleaning products - such as benzyl alcohol, alcohol wipes or liquid detergents.